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Opracowano: 2019-03

Instruction for use
Read the instruction carefully before use.

Endostar Gutta Cleaner

Gutta-percha orange solvent

FOR PROFESSIONAL DENTAL USE ONLY

Intended use

Product for dissolving gutta-percha filling during endodontic retreatment.

It is especially recommended:

- in case of difficult gutta-percha removal from curved root canals,
 - for very well condensed and hardened gutta-percha and sealer removal.
- Helpful while cleaning of dental instruments from gutta-percha and sealer remnants.

Technical and medical data

Composition: sweet orange oil 100%.

Primary package: amber glass bottle with dropper, contains 10 ml of the liquid product.

Mechanism of action

Endostar Gutta Cleaner is a solvent. It dissolves gutta-percha and root canal sealer (particularly zinc oxide-eugenol based sealer), making them easier to remove from the root canal while endodontic retreatment.

Directions for use

Gutta-percha and root canal sealer removal

- using rotary instruments (e.g. Gates-Glidden drill) remove filling material from tooth chamber, creating reservoir for the solvent,
- place 1-2 drops of solvent into established space,
- after approximately 2 minutes progressively remove softened gutta-percha with endodontic hand file e.g. K-file, adequate to root canal diameter (e.g. ISO 15 or 20) or H-type file,
- if necessary use next dose of the product, place it in partially cleaned root canal,
- during removal procedure irrigate root canals with 2-3% sodium hypochlorite solution,
- stop using solvent when the endodontic file reaches the distance 3-4 mm from root apex,
- irrigate root canals accurately (e.g. with 2-3% sodium hypochlorite solution),
- remove remaining filling material using endodontic files and irrigate the root canal thoroughly with 2-3% sodium hypochlorite solution,
- follow endodontic retreatment protocol afterwards.

Cleaning of dental instruments from gutta-percha and sealer remnants

- clean instrument with gauze soaked with solvent,
- rinse instrument with warm running water.

Warnings and precautions

- ☞ Do not use the product in the case of hypersensitivity to sweet orange oil.
- ☞ Do not use the solvent to remove gutta-percha extruded beyond the root apex.
- ☞ Do not let the solvent to infiltrate beyond the root apex.
- ☞ Use rubber dam during application of the product. It protects soft tissues of oral cavity from contamination with the product.

Causes skin irritation.

May cause an allergic skin reaction.

May be fatal if swallowed and enters airways.

Wear protective gloves/protective clothing/eye protection/face protection.

IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

IF INHALED: If experiencing respiratory symptoms: Remove person to fresh air and keep comfortable for breathing.

Flammable liquid and vapour.

Keep away from heat, hot surfaces, sparks, open flames and other ignition sources.

No smoking.

Very toxic to aquatic life with long lasting effects.

Dispose contents/container in accordance with requirements of utilization.

Place and storage conditions

Keep container tightly closed. **Keep out of the reach of children.**

Store in temperature below 25 °C.

Keep away from sunlight.

The lot number is given on the outer and primary package.

The expiry date is given on the outer and primary package. Do not use the product after expiry date. After opening the expiry date does not change provided that the container is tightly closed after each use.

Explanation:

☞ - pay special attention to the texts marked with this symbol

Prepared: 2019-03